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Copyright © 2013, GREINER GmbH
All parts of this documentation are protected by copyright. Any utilisation or modification, which falls outside the close restrictions of the Copyright Act is not permitted without the agreement of GREINER GmbH and will be liable to prosecution. This applies in particular to copying, translating, microfilming and to storing and processing in electronic systems.

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
Händlerstempel / dealer stamp

Contents

1	Equipment safety	4
1.1	Declaration of conformity	4
1.2	Storage and transport information	4
1.3	Information for setting-up and use	4
1.4	Requirements for use	5
1.5	Warning and information labels	5
1.6	Instruction	5
1.7	Safe operation	5
1.8	Meaning of symbols	6
2	Product identification	6
2.1	Model variants	6
3	Description and operation	7
3.1	General	7
3.2	Operating the components	7
3.3	Electrically adjustable models	8
3.4	Charging the battery	9
3.5	Changing the seat cushion	9
3.6	Operating the leg support	9
3.7	Length adjustable footrest	9
3.8	Locking the castors	10
3.9	Pivoting foot rest	10
4	Maintenance.....	10
4.1	General	10
4.2	Mechanical components	10
4.3	Testing electrical safety	10
5	Technical data	11
6	Care.....	12
7	Additional information for models with electrical adjusting elements	12
8	Accessories	13
9	Spare parts overview	14
10	Disposal.....	14
11	Maintenance checklist	15

1 Equipment safety

The GREINER carryLine **mobile** **lounger** has been designed and tested in compliance with national and international specifications and is manufactured in accordance with the applicable directives. This ensures that the required level of safety is provided by the manufacturer.

Please read these instructions for use carefully and make yourself familiar with the safety information before using the GREINER carryLine **mobile** **lounger**.

Safety information is indicated by the adjacent symbol and demands your particular attention!

Please also observe the instructions for use for any other equipment and accessories. You can obtain further information from our Service Department or from authorised representatives.

1.1 Declaration of conformity

The GREINER carryLine **mobile** **lounger** described in this manual has been designed and manufactured in accordance with EC Directive 93/42 EEC. It is a medical product of Risk Class 1 as defined by the Medical Product Act (MPG).

The Declaration of Conformity shall cease to be valid if changes are carried out to the product by the customer or third parties. Examples of such changes are modifications of all kinds, use of third-party accessories and the removal of warning and information labels.

1.2 Storage and transport information

The carryLine **mobile** **lounger** may only be transported in an upright position and carefully secured.

The following environmental conditions must be maintained for the storage and transportation of models with electrical adjustment devices:

Temperature	+10°C to +50°C
Humidity	20% to 90% at 30°C
Air pressure	700 to 1,060 hPa

1.3 Information for setting-up and use

The carryLine **mobile** **lounger** is intended exclusively for use inside buildings under normal environmental conditions.

The following environmental conditions must be maintained when using models with electrical adjustment devices:

Temperature	+5°C to +40°C
Humidity	20% to 90% at 30°C
Air pressure	700 to 1,060 hPa

Intended use:

- in examination and treatment rooms
- in homes for the aged and nursing homes
- care facilities
- hospitals
- transportation of patients
- Intended users: Patients under survey of trained and instructed personnel of the concerned department. In addition: doctor, assistant, sister, curator, etc. .
- Installation and commissioning: by the user's technical personnel or our contract dealers.

NON-intended use:

- In the vicinity of high-frequency surgical units or defibrillators where these do not comply with the requirements of EN 60601-1-2:2001
- In an environment where the occurrence of inflammable or explosive gases or vapours (e.g. anaesthetics) is to be expected. Applies to electrically operated models.
- Work chair
- Operating chair / table
- Alternative to a bed
- Emergency transport

1.4 Requirements for use

The following points must be observed in order to guarantee safe operation:

- Make yourself familiar with the content of these operating instructions before using the GREINER carryLine **mobile** lounge. In particular, pay attention to the safety information relating to the individual areas of risk.
- The lounge must be disinfected before being used for the first time and thereafter as required.
- All clamping levers and clamping screws, including those on accessories, must be tightened before use.
- Access to the mains plug must be guaranteed at all times.

The lounge must not be cleaned in a washing machine or with a high-pressure cleaner. Spraying down with a water jet is also not allowed.

1.5 Warning and information labels

Warning and information labels can be found in the following positions:

- nameplate with CE mark and serial number and all technical information at the right-hand side of the metal frame below the seat.
- Nameplates with appropriate data on all electrical components.

If you should establish that one of these labels is missing from your equipment or has become illegible, please contact your dealer. He will provide you with a replacement.

1.6 Instruction

A medical product may only be put into operation and applied if the manufacturer, or a person authorised by the manufacturer, has first subjected the medical product to a *functional check* at the place of use. Furthermore, the persons authorised by the user must be instructed by the manufacturer, or a person authorised by the manufacturer, in the handling, use and operation of the medical product based on the instructions for use and the attached safety information and commissioning instructions. The instruction must be documented. You will find an appropriate confirmation in the appendix.

1.7 Safe operation

If the GREINER carryLine **mobile lounge is equipped with electric adjusting drives and extremely sensitive measuring instruments are used at the same time, electromagnetic distortion of the measurement results while the mains plug is connected and/or the adjustment motors are running cannot be ruled out!**

Equipotential bonding must be used if other electrical medical equipment is to be used with the chair (see 4.3). Equipotential bonding is a preventative measure for protecting patient, user and third parties against contact voltages, which can occur due to voltage potentials between accessible conducting parts and mobile units.

The GREINER carryLine **mobile** lounge must always be maintained in perfect and fully functional condition. You must not use the lounge if a defect occurs, which you are unable to rectify. Mark the equipment as being non-operational and inform our Service Department.

Mechanical stressing of the mains cable must be avoided. Regular visual inspection of the mains cable at short intervals is required, particularly after any mechanical stressing. The mains cable is not resistant to being driven over! Damaged cables may only be replaced by an electrical specialist. Please note the additional information for models with electrical adjustment devices and batteries in Chapter 7.



Always move the lounge to the normal position to get in and out: seat set to the lowest height, backrest upright, legrest retracted. To transport the patient the lounge can be declined slightly, a transport of the patient in the flat position is not allowed. Disconnect the mains cable first and do not push the lounge over the mains cable.

Do not sit on the leg rest or backrest. There is a risk of tipping if the chair is additionally or incorrectly loaded!

This must be expressly pointed out to the patients as well as to the care personnel and users!

1.8 Meaning of symbols

-  CE Declaration of Conformity
-  Caution! Safety information
-  Read the operating instructions
-  For indoor use
-  Protection class II
-  Application part Type B
-  Alternating current
-  Do not dispose of with residual waste!
-  Manufacturer
-  Date of manufacture
-  Attention, please refer to accompanying documents
- SN** Serial number
- REF** Item number
-  Equipotential bonding
-  Connection for equipotential bonding
-  No sitting
-  Maximum permissible loading

2 Product identification

If damage or a defect should occur with your GREINER carryLine **mobile lounger** or if you should have any other reason for complaint, please have a look at the inside of the right-hand side of the matl frame below the seat. Here you will find the product sticker shown below, which contains the following information:

- Date of manufacture
- Item number and description
- Serial number
- AB number
- Customer number

		GmbH Riedbachstr. 5 D-74385 Pleidelsheim
REF:	4870200	  
Polster / Cushion:	A15	
SN:	20094567800050001	 200kg
	 01.02.09	
Kunde / Customer:	020046 Klinik Musterhausen	
U in:	240V ~, 50/60Hz	I in: max. 0,25A
	24 V 	
Op:	5%, max. 3min/57 min.	
Symbol:	   	Q-Control: 1033
Code:		

Please give this information to your dealer or to our Service Department in case of any complaint. The more accurate and better the information is, the better and more objective will be the help provided.

2.1 Model variants

The GREINER carryLine **mobile lounger** is available in different configurations. These operating instructions apply to all those models listed in the following matrix.

	Adjustment		Foot rest	
				
Modell	manual	electrical	pivoting	length adjustable
4870200	•			
4870210	•		•	
4870300		•		
4880310		•	•	
4870400	•			•
4870500		•		•

3 Description and operation

3.1 General

All clamps, including those of the accessories, must be checked for tightness before each use.

The equipment must only be operated by trained personnel.



The GREINER carryLine mobile lounger must always be moved to the basic position for getting in and out – base unit at the lowest position, upper section upright. Make sure that the patient's hands are on the armrests when adjusting the lounger.

The armrests (2) can be folded out of the way when the locking bolt (3) is pulled down. When folding the armrests (2) forward again, ensure that the locking bolt (3) re-engages correctly.

3.2 Operating the components

The back and leg supports are adjusted synchronously with one another. To do this, push the release lever (1) back and keeping it pressed bring the back section (4) into the required position and then let go of the release lever (1). Repeat the procedure in the reverse order to move the back (4) into an upright position once more. Stepless adjustment of the back is achieved by a lockable gas strut.

With electrically adjustable models, press the appropriate button on the manual controller (1b).

3.3 Electrically adjustable models

Caution: Models with electrical adjustment must not be used in areas where there is a risk of explosion!

All electrically adjustable models described here are operated with an intermittent duty of 3 min/57 min. i.e. run for max. 3 minutes at rated load followed by a pause of 57 minutes. Otherwise, functional failure may occur!

All electrically adjustable positions (5) can be blocked by the enclosed key (6) and the blocking buttons (7) to avoid an inadvertent activation of the functions.

Make sure that no one is in the danger area when operating the GREINER carryLine *mobile* lounger


à Danger of getting caught! Please also advise the patient.

The GREINER carryLine *mobile* lounger is not intended for use by small children or frail persons without supervision. Keep the manual controller out of reach when there are no trained person in the vicinity of the lounger.


3.4 Charging the battery

Before the GREINER carryLine mobile lounger is used for the first time, the battery must be charged for 24 hours in order to achieve correct operation and an extended life. To charge the battery, the power supply unit (1) has to be plugged in the socket (2) at the right side part. You will find further information on the battery in Chapter 7 on Page 12.

3.5 Changing the seat cushion


The seat cushion (1) can be removed for maintenance or cleaning purposes. To do this, raise the front of the seat cushion (1) until the fixing clip (3) disengages from the bar (4), and then remove the seat cushion by pulling forwards.

When replacing, make sure that the metal tongue (5) is positioned underneath the mobilepiece (6), and the fixing clip audibly re-engages on the bar (4). The fixing bridge (2) makes positioning easier.

3.6 Operating the leg support

The leg support (1) can be adjusted independently of the synchronous movement of the GREINER carryLine mobile lounger. To do this, release the retaining clip (3) from the frame (4) with a brief tug, raise the leg support (1) and release it in the required position. The latching mechanism (2) stays at the required height.


To lower the leg support (1) once more, raise it slightly and then lower it smoothly and engage the retaining clip (3) with the frame (4).

3.7 Length adjustable footrest

In models 4870400 & 4870500 the legrest moves together with the back of the lounger. The extension position of the footrest (1) can be changed.

To do this, lift the footrest frame (1) on both sides from


its slots in the guide tubes (2) and insert it in a different position.

Do not step on the footrest (1) when sitting down or getting up. Danger of tipping! The footrest (1) is designed only to support the patient's feet.

3.8 Locking the castors

The GREINER carryLine **mobile lounger** is equipped with two castors. The castors can be blocked by stepping on the locking device (1).

Never leave the GREINER carryLine mobile lounger on ramps or sloping surfaces.


The castors must always be locked while the lounger is not in use as the chair may suddenly start rolling!

During patient transport make sure there is always sufficient space between the user and walls. Risk of bruise injury!

3.9 Pivoting foot rest

The foot rest can be pivoted underneath the chassis frame. Therefore push the step plate (1) with foot or hand. Ensure that the stop bolt (2) engages properly in the spring clamps (3). The foot rest is lowering during stepping on it so the danger of tipping and stumbling can be avoided.

4 Maintenance

4.1 General

As an operator as defined by the Medical Product Act (MPG) you are obliged to protect users, patients and third parties against hazards that can arise from medical products by means of regular inspections and checks. We therefore recommend that repair and maintenance work as well as safety checks in accordance with §4 of the ordinance concerning the setting up, operation and application of medical devices (MPBetreibV) be carried out and recorded at regular in-

tervals. It is the responsibility of the operator to carry out inspections and to define the inspection intervals.


4.2 Mechanical components

Check all screws, hand wheels, clamp levers, castors etc. for completeness and tightness. Missing fixings must be replaced and loose screws must be tightened. Check the completeness of all axles and whether they are all secured against dropping out by means of two circlips. Treat all moving linkages, bearings and motor piston rods with suitable lubricants, e.g. silicone spray.

4.3 Testing electrical safety

Electrical safety must be checked in accordance with the specifications of the employer's liability insurance association for your area (e.g. BGV A3). Recurring measurements have to be done according to the enclosed test protocol (electrical safety) To carry out a visual inspection, remove the seat cushion as described in Chapter 3.5.

Check all cables for breakage and chafing and for crushing. If a cable is defective, it must be replaced before the inspection can be continued. Check the housing of electrical components such as the battery controller (1), motor (2), charger (3) and manual controller (4) for cracks and breaks. Make sure that the equipotential bonding pin (5) is tight and not bent or otherwise damaged.

Caution: Work on the electrical system and/or the safety checks must only be carried out by an authorised customer service engineer or by a trained specialist.

You will find a checklist with the minimum requirements for regular maintenance in the Appendix.

5 Technical data

Mechanical data:

Seat width:.....	486mm
Seat depth:.....	ca. 500mm
Overall height (without accessories):.....	1295mm
Seat height:.....	503mm
Overall width (without accessories):.....	696mm
Overall height (in lying position):.....	928mm
Lift height ca.:.....	190mm
Overall depth (in lying position):.....	1667mm
Seat angle with tilt mechanism.....	ca. 10°
Maximum permissible loading:.....	200kg

Adjustment motor data:

Rated voltage.....	24V DC
Current consumption at rated load.....	max. 5A
Thrust:.....	3000N
Degree of protection.....	IP51

Battery controller data:

Rated voltage.....	24V DC
Capacity.....	2.7Ah
Degree of protection.....	IPX0

For indoor use:..... 

Protection class II..... 

Application part Type B..... 

Mains connection.....230 V AC 50/60 Hz




Operating mode¹⁾ Intermittent duty 3min./57 min.

1) Operating mode = Intermittent duty 3 min./57 min. i.e. run for max. 3 minutes at rated load followed by a pause of 57 minutes. Otherwise, functional failure may occur!

Risk of tipping:

The maximum load of 200kg must be distributed evenly over the lounge surface. Do not sit on the leg support, arm-rest or backrest!

6 Care

The GREINER carryLine *mobile* lounge must be cleaned and disinfected in accordance with the general basic principles before using for the first time and after every use.

The plastic parts used are made from ABS or integral foam and must not be cleaned with petrol, benzene, turpentine or other solvents. Lukewarm soap solutions or commercially available plastic cleaners are most suitable for cleaning plastic parts. Follow the instructions for use for the respective cleaner.

On no account must cleaning agents, which contain abrasive particles, be used for cleaning steel parts.

The cleaning, maintenance and care of electrical components is described in Chapter 7.

The cushions materials must be protected from the effects of excessive heat and they need to be cleaned at regular intervals. When cleaning the cushions, care must be taken to ensure that they do not become saturated. Slight soiling can be wiped away with a damp microfiber cloth. Stain caused by oil, fat or sweat must be removed immediately. Therefore use a warm and mild soap sud and a microfiber cloth or a soft brush. Do not use solvents, chlorides, abrasives, chemical cleaning agents or wax polishes. Avoid any contact with cleaning agents which contain oil, fat or alcohol. The upholstery material is not resistant to chemical and dry cleaning.

Caution: The GREINER carryLine *mobile* lounge must not be cleaned in a washing machine or with a high-pressure cleaner.

When disinfecting, do not allow the disinfectant to soak into the cushion for more than the specified time. Disinfectants are aggressive and can cause a change in the surface. Use disinfectants only as recommended by the manufacturer and in the appropriate concentration.

In order that the operation of moving parts is not detrimentally affected, these must be cleaned regularly, axles and linkages must be lubricated with commercially available oil or grease, and fixing screws must be checked for tightness.

Caution: Maintenance work, which extends beyond the described scope, must only be carried out by an authorised customer service engineer or by a trained specialist.

7 Additional information for models with electrical adjusting elements

Before installation, de-installation or troubleshooting:

- Pull the mains plug.
- Remove all loads from the adjusting drives. System connection: the individual components must be connected before the control box is connected to the mains.
- Ensure that the voltage stated on the control box is correct before connecting the system to the mains.
- When connecting the control box, ensure that the cable is not trapped and that it is strain-relieved, and also that it is not damaged by sharp-edged parts during the various directional movements.

During operation:

- Switch off the power supply immediately if unusual noises or odours come from the control unit.
- Ensure that the cables are not damaged.
- Pull out the mains plug before moving any mobile parts.
- The products may only be used in an environment matching their IP protection class.
- The duty cycle stated on the control box (see chap. 5 page 11) must always be maintained. If it is exceeded, overheating of the control box may occur. The electrical components are not designed for use in areas where there is a risk of explosion!

Care and maintenance:

- The materials are resistant to the majority of cleaning agents and disinfectants as used in hospitals and in domestic nursing.
- Cleaning with a high-pressure cleaner or steam cleaner is not permitted, as it is difficult to maintain the necessary distance of 300 mm from the electric components.
- To prevent degreasing of the piston rods, the adjusting drives should be returned to the smallest stroke setting before cleaning.

Disinfectants and cleaning agents:

- The agents must have a pH value that is neither too alkaline nor too acid (6 – 8).
- They must not contain any caustic constituents. They must be free of additives, which can attack the surface or can dissolve plastic parts.
- They must not attack lubricating greases.

Battery systems:

- Batteries must be replaced after no later than 4 years, and possibly earlier, depending on use. Frequent and powerful discharging will reduce the service life of the battery.
- For an optimum service life, it must be connected as often as possible to the mains voltage.
- Batteries must be recharged at least every 3 months, otherwise they will be damaged by self-discharging.
- The operation of the battery must be checked once a year
- Batteries may only be replaced by mechanically and electrically compatible types. A replacement can be obtained from Greiner or from your dealer.
- When fitting replacement batteries, ensure that they are connected correctly. Compare with the drawing in the battery compartment and check that no connections are loose.
- **CAUTION:** In the factory, the battery compartment is hermetically separated from the electronics area. When batteries are replaced, this separation must not be damaged or modified, as in this way it is possible for battery gases to penetrate the electronics area, which in turn would present a risk of explosion!
- **Disposal:** Lead-acid batteries can be returned to your dealer or disposed of in the same way as car batteries.
- **CAUTION:** The battery compartment is provided with ventilation. This must not be blocked or covered, as an explosion may result due to over-pressure.
- If the battery housing has been mechanically overloaded (dropped on the floor, collision/compression in an application or heavy blow), it must be sent to an authorised service workshop for checking the hermetic separation between battery compartment and electronics compartment.

8 Accessories

Cushions:

Additionally a pair of supporting cushions for the back section (2), and the tights (1), and different neck supports (3) are available.

Storage basket:

The storage basket (4) is used for storing the patient's personal belongings for example.

Seatbelt

Caution: The seatbelt (5) is used for securing the patient during transportation. Under no circumstances may it be used for securing unattended persons.

Patient table

The patient table (6) is clamped to the armrests using the clamps fitted to both sides.


Risk of scalding! Do not serve hot meals or drinks on the table.

Intravenous drip stand:

The intravenous drip stand (7) is adjustable in height.

Rails:

Standard rails (8) for fixing specific equipment are also available.

9 Spare parts overview




Pos:	Art.Nr.:	Bezeichnung:
01	4879273	tension spring
02	4869201	pneumatic spring
03	4879275	bowden wire
04	4439203	release button
05	4879913	pushing handle
06	4877803	backrest cushion
07	4877703	seat cushion
08	4877403	leg support cushion
08b	4877404	leg support cushion, height adjustable*
09	4879215	wheel Ø300 PU-tires
10	4879214	castor Ø125mm
11	4480115	latching mechanism
12a	4879217	armrest integral
13	4879229	armrest locking bolt
14	6409220	adjustment drive
15	4469231	battery box
16	4469239	battery charger CH01
17	4469230	manual control panel
18	4469236	Key for manual control panel

*only in models with height adjustable footrest 4870400/4870500

10 Disposal

When the goods supplied are electrical or electronic units as defined by the stipulations of the act governing the Sale, Return and Environmentally Sound Disposal of Electrical and Electronic Equipment (ElektroG) of 16.03.2005 and DIRECTIVE 2002/96/EC OF THE EUROPEAN PARLIAMENT AND THE COUNCIL of 27th January 2003 concerning Waste in Electrical and Electronic Equipment (WEEE) (e.g. examination loungers with electric motor adjustment), GREINER will take back the goods supplied at the end of their useful life at the customer's expense and dispose of them in an orderly manner in accordance with the statutory requirements. The customer is obliged to make the equipment available to GREINER in an orderly manner. The equipment is identified by means of the serial number and the GREINER item sticker on the unit. It is incumbent upon the customer to show that the unit is a GREINER product.

Please refer to our General Conditions of Business at www.greiner-gmbh.de for more information

11 Maintenance checklist

- All mechanical fixings tight?
- All screws tight, all axles present and fitted with two lock washers?
- Equipotential bonding pin present and undamaged?
- Name plate present and legible?
- Hand wheels, clamping levers and locking knobs complete and operational including those of the accessories?
- All Bowden cables intact, wires not split? Oil Bowden cables sparingly.
- Oil leaking from gas struts? If yes, have them replaced!
- All welds intact, not cracked or damaged?
- All castors secured with screws, rotate and swivel without problems, holding brake working?
- Upholstered parts all undamaged, no burst seams, no tears or cuts?
- Mains cable and all other cables (charger, manual controller, battery controller, motor) undamaged, no breaks, kinks or chafed areas?
- All electrical component housings undamaged? No breakages or other damage?
- Electrical safety checked and documented?
- Check motor piston rods. If necessary, grease with a thin coating of moisture-free Vaseline.
- No bent, damaged or sharp-edged parts.
- Seatbelt (accessories): Belt not damaged or frayed, all seams OK, seatbelt lock working?

Serial number of chair: _____ Model No.: _____

Maintenance work carried out by: _____ Date: _____

This checklist is intended as an aid in meeting the minimum requirements of regular maintenance. It therefore does not constitute a conclusive, comprehensive overview. Type, scope and repetition of the maintenance work to be carried out depend on the application and must be assessed by the operator.

Confirmation of instruction

Models 487-488 carryLine

Riedbachstrasse 5
D-74385 Pleidelsheim
07144 / 8112- Tel: -0 Fax: -99
info@greiner-gmbh.de

In accordance with §5 of the **MPBetreibV** (German ordinance concerning the setting up, operation and application of medical devices), the medical product which is subject of this instruction may only be put into operation and applied if the manufacturer, or a person authorised by the manufacturer, has first subjected the medical product to a functional check at the place of use. Furthermore, the persons authorised by the user must be instructed by the manufacturer, or a person authorised by the manufacturer, in the handling, use and operation of the medical product based on the instructions for use and the attached safety information and commissioning instructions. The instruction must be documented. When the persons authorised by the user have been familiarised with the medical product and the operating manual, the trainer must once again clearly emphasise the following particularly relevant points:

- The lounge must not be loaded with more than 200kg.
- In spite of EMC testing (electromagnetic compatibility), and this has occurred on previous occasions, measuring equipment connected to the patient may be affected when the electric motors of the lounge are operated while measurements are being taken.
- Use equipotential bonding when measurements or treatment are carried out on the patient with other electrical medical products.
- Always move the lounge to the normal position to get in and out: seat set to the lowest height, backrest upright, legrest retracted. To transport the patient the lounge can be declined slightly, a transport of the patient in the flat position is not allowed. Disconnect the mains cable first and do not push the lounge over the mains cable
- **carryLine cross**: Only cross thresholds in reverse. If a patient is seated in the lounge, it will be impossible to cross the threshold. If an attempt is made to do this by building up momentum, this could cause damage to the lounge and, in the worst case, it could tip.
- **carryLine hvh & carryLine mobile**: Do not step on the length adjustable footrest when sitting down or getting up. Danger of tipping! It is designed only to support the patient's feet.
- Do not sit on the leg rest or back section as otherwise the lounge could tip.
- Always move the lounge to the normal position when patients get in or out.
- The brakes must be set on the castors when the lounge is parked with or without the patient. → Otherwise, the lounge could unintentionally roll away.
- Caution: The **carryLine** lounge must not be used in areas where there is a risk of explosion. The controller and motors are not explosion proof. In addition, in battery mode, it must be noted that the battery releases gases which may possibly form an explosive mixture with other gases.
- The electrical components are only proof against splash water. So the lounge must not be cleaned in a washing machine or with a high-pressure cleaner.
- Regularly check axles, bearings and articulated joints for abrasion and wear. Replace worn parts. All circlips must be fitted at all times.

Instruction carried out by:	The following were given instruction:
Date:	1.)
Name, address, company stamp:	2.)
	3.)
	4.)
	5.)