Endosgel® Chlorhexidindigluconat 0,05% DE

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Sans latex.

Patient Information Leaflet**Endosgel® Chlorhexidine Digluconate 0,05% EN**

Sterile lubricant with physiological conductivity for use in electro-resection and electro-coagulation and for intermittent catheterization.

Contains: Chlorhexidine Digluconate 0,05%, Methyl Hydroxybenzoate, Propyl Hydroxybenzoate, Hydroxyethylcellulose.

Properties: Endosgel® is a sterile lubricant with an electrical conductivity in the physiological range. Endosgel® adheres well to the mucous membrane, has excellent lubricant properties and maintains a clear view through the optics of instruments. Endosgel® has auxiliary disinfectant properties.

Indications: Electro-resection and electro-coagulation, transurethral endoscopic surgery, transvaginal endoscopic surgery, stricturotomy, resection of the neck of the bladder, prostate resection, tumor coagulation, resection of the bladder, intermittent catheterization (6 or 11ml), endourological surgery, rectal surgery. **5-10 minutes should elapse for the gel to reach its full effect.**

Contraindications: Must not be used in patients who are hypersensitive/allergic to alkyl hydroxybenzoates (parabens), chlorhexidine or any of the other constituents.

Side effects: Hypersensitivity may result in skin symptoms; in very rare cases, contact dermatitis or irritation of the bladder may occur. In very rare cases (less than 0.01%) allergic reactions to chlorhexidine digluconate can occur, including anaphylactic shock.

Interactions with other medicaments and other forms of interaction: Until now none known.

Special precautions for use: In view of the chlorhexidine content of the product, prolonged use of Endosgel® is inappropriate unless a doctor has been consulted.

Dosage instructions: Unless otherwise prescribed, expel 6/11 or 20ml Endosgel® from the single-use syringe into contact with the mucous membrane (instill). Further doses may be applied if necessary.

Urethral use: Perform the usual cleaning and disinfection of the surrounding skin and/or mucous membrane (e.g. the external genitalia and the orificium urethrae externum) and then instill 6/11 or 20ml Endosgel® into the urethra in a slow and measured dose. The dose depends on the patient's anatomy and/or the type of intervention to follow: Single injection of 6ml Endosgel® is suitable for use in women and children; 11ml is suitable for men, and 20ml is recommended for electro-resection and electro-coagulation.

Other uses: 6/11 or 20ml Endosgel® can be evenly applied to the mucous membrane or to the instrument to be lubricated (e.g. catheter or endoscope).

Advices and information concerning the stability of the medical device: Endosgel® must not be used after the expiry date. Discard remaining gel and do not use the empty syringe for any other purpose to avoid the risk of transmission of pathogens.

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Medeco B.V.
Germania Pharmazetika GmbH
Melisana AG
AD CARE AB
FARCO-PHARMA GmbH